

برگشت (/Panel/Labs/Answers/47998)

چاپ

دوره پنجاه و سوم زمستان 1404

S2 Results

نمونه مجهول ToRCH

کد نمونه 1173	EQAP-18421	آزمایشگاه آزمایشگاه آرامش
---------------	------------	---------------------------

Anti-CMV (IgG)

Anti-HSV1 (Anti-Herpes Type1 IgG)

Anti-Rubella (IgG)

Anti-Toxoplasma (IgG)

Anti-Rubella (IgG)

جدول شماره 1. میزان مشارکت در این دوره-آزمایش

جدول شماره 2. توزیع فراوانی نتایج گزارش شده توسط شرکت کنندگان براساس روش

تذکر. نتیجه مورد قبول برای این آزمایش پاسخ + بوده است. نتیجه ارائه شده توسط آزمایشگاه شما با علامت

← مشخص شده است.

تعداد ثبت نام کنندگان	تعداد شرکت کنندگان	کل درصد مشارکت	تعداد پاسخهای دریافتی	درصد پاسخدهی
191	168	% 87.96	161	% 95.83

Method	Count				Percent		
	Total	Positive	Indeterminate	Negative	Positive	Indeterminate	Negative
CL	48	47	0	1	%97.92	%0.00	%2.08
ECL	24	24	0	0	%100.00	%0.00	%0.00
ELFA	2	2	0	0	%100.00	%0.00	%0.00
ELISA	79	77	0	2	%97.47	%0.00	%2.53
Immunochromatography	0	0	0	0	%0.00	%0.00	%0.00
Other	6	6	0	0	%100.00	%0.00	%0.00

جدول شماره 3. توزیع فراوانی نتایج کیفی گزارش شده به تفصیل بر اساس کیت و روش

Method	Kit	Result	Count	Your Answer
CL	Abbott	Positive	11	
	Architect	Positive	2	
	Cobas	Positive	1	
	Dia Sorin	Negative	1	
		Positive	20	
	Liaison	Positive	2	
	Maglumi	Not Reported	1	
		Positive	4	
	Ortho	Positive	1	
	Pishtaz Teb	Positive	1	
	Roche	Positive	1	
Siemens	Positive	4		
ECL	Not Reported	Not Reported	1	
	Abbott	Positive	1	
	Cobas	Positive	2	
	Liaison	Positive	1	
	Maglumi	Positive	1	
	Roche	Positive	18	
	Siemens	Positive	1	
ELFA	Biomerieux	Positive	1	
	Chorus	Positive	1	
ELISA	Not Reported	Positive	2	
	Chorus	Positive	5	
	Diesse-Tec	Positive	2	
	Euroimmun	Positive	3	
	Golden	Positive	1	
	Pishtaz Teb	Not Reported	1	
		Negative	2	
		Positive	60	←
Vircell	Positive	4		

Method	Kit	Result	Count	Your Answer
Other	Chorus	Positive	3	
	Diesse-Tec	Positive	2	
	Maglumi	Positive	1	

Copyright © 2026 IACLD. All rights reserved.